

Num	Product Name	Pore Size	Surface area (m ² /g)	End-capped	C%	Optimum pH	pH Range
1	Unisol C18	100Å	410	Single	18	2.0-6.0	1.5-9.0
2	Unisol C8	100Å	410	None	13	2.0-6.0	1.5-9.0
3	Unisol C18(2)	110Å	340	Single	14	2.0-6.0	1.5-9.0
4	Unisol Amide	100Å	410	None	8	2.0-6.0	2.0-8.0
5	Venusil AQ C18	100Å	380	Single	18	2.0-6.0	1.5-9.0
6	Durashell C18-AM	100Å	380	Double	14	2.0-10	1.5-12
7	Durashell C18	100Å	380	Double	21	2.0-10	1.5-12
8	Durashell C8	100Å	380	Single	14	2.0-10	1.5-12
9	Durashell RP	150Å	200	Double	15	2.0-10	1.5-12
10	Durashell C18(L)	150Å	200	Double	14	2.0-10	1.5-12
11	Durashell NH ₂	100Å	380	None	9	2.0-8.0	2.0-10
12	Venusil ASB C18	150Å	200	None	12	1.0-4.0	0.8-7.5
13	Venusil ASB C8	150Å	200	None	7	1.0-4.0	1.0-7.5
14	Venusil ASB C1	150Å	200	None	3	1.0-4.0	1.0-7.5
15	Venusil ASB Phenyl	150Å	200	None	6	1.0-4.0	1.0-7.5
16	Innoval C18	100Å	220	Double	14	2.0-6.0	1.5-9.0
17	Innoval ODS-2	100Å	220	Double	14	2.0-8.0	1.5-9.0
18	Innoval AQ C18	100Å	220	Single	15	2.0-6.0	1.5-9.0
19	Innoval Silica	100Å	220	None	-	2.0-5.0	1.0-7.0
20	Venusil XBP C18	100Å	380	Double	22	2.0-6.0	1.5-10.0
21	Venusil XBP C18(2)	100Å	380	Double	18	2.0-6.0	1.5-9.0
22	Venusil XBP C18(A)	120Å	340	Double	20	2.0-6.0	1.5-9.0
23	Venusil XBP C18(L)	150Å	200	Double	15	2.0-6.0	1.5-9.0
24	Venusil XBP C18(T)	300Å	80	Double	8	2.0-6.0	1.5-9.0
25	Venusil XBP C8	100Å	380	Single	14	2.0-6.0	1.5-8.5
26	Venusil XBP C8(L)	150Å	200	Single	8	2.0-6.0	1.5-8.5
27	Venusil XBP C4	100Å	380	None	7	2.0-6.0	1.5-8.5
28	Venusil XBP C1	100Å	380	None	3	2.0-6.0	1.5-8.5
29	Venusil XBP NH ₂	100Å	380	None	6	2.0-6.0	2.0-8.5
30	Venusil XBP Phenyl	100Å	380	None	12	2.0-6.0	1.5-8.5
31	Venusil XBP CN	100Å	380	single	5	2.0-6.0	2.0-9.0
32	Venusil XBP Silica	100Å	380	None	-	2.0-5.0	1.0-7.0
33	Venusil XBP Silica(L)	150Å	200	None	-	2.0-5.0	1.0-7.0
34	Venusil SCX	150Å	200	None	6	2.0-6.0	2.0-8.0
35	Venusil SAX	100Å	380	None	6	2.0-6.0	2.0-8.0
36	Venusil XBP COOH	100Å	380	None	3	2.0-6.0	2.0-8.0
37	Venusil Diol	100Å	380	None	-	2.0-6.0	2.0-8.0
38	Venusil AA	100Å	380	Double	22	2.0-6.0	1.5-9.0
39	Venusil PAH	200Å	200	None	17	2.0-6.0	1.5-9.0
40	Venusil PFP	120Å	320	None	-	2.0-6.0	2.0-8.0
41	Venusil HLP C18	110Å	340	None	15	2.0-6.0	2.0-8.0

42	Venusil XBP Polar- Phenyl	120Å	340	Single	11	2.0-6.0	1.5-9.0
43	Promosil C18	100Å	320	Double	18	2.0-6.0	1.5-9.0
44	Promosil C8	100Å	320	Double	10	2.0-6.0	1.5-9.0
45	Promosil Silica	100Å	320	None	-	2.0-5.0	1.0-6.0
46	Promosil NH2	100Å	320	None	4	2.0-5.0	2.0-8.0
47	Venusil PrepG C18	120Å	340	Double	17	2.0-6.0	1.5-9.0
48	Optimix C18/C8	100Å	380	None	-	2.0-6.0	1.5-9.0
49	Optimix C18/SCX	300Å	80	None	-	2.0-6.0	1.5-9.0
50	Optimix C18/Amide	100Å	380	None	-	2.0-6.0	1.5-9.0