

Vortex Mixer | AVM2001



제품 규격 및 사양

Model	AVM2001
Type of movement	Orbital
Speed Range	0 to 3,300 rpm
Operating Modes	Touch or continuous
Orbit diameter	Circular Orbit, 4.0mm Dia
Speed Control	Stepless
Application Container / Number	Hand holding; any micro tube or conical tube Fixed holding (Micro tube Platform Head) : 0.5/1.5/2.0ml micro tube x 12 each
Construction Material	Powder Coated Aluminum
Driving Motor	Brushless shade motor
Permissible Temp. & Humidity	5℃~40℃, <80%
Max. Load & Power	0.5Kg, 30W
Power Supply	AC 120V, 60Hz or AC 230V, 50/60Hz



AVM2001RP



AVM2001PC



AVM2001RA

기본구성품

Ordering No.	Description
AVM2001RP	3inch Rubber Platform Head
AVM2001PC	Pop-off Cup Head
AVM2001RA	Remove Applicator

* 제품의 규격 및 사양은 사전공지 없이 변경될 수 있습니다. 발주전 문의 부탁드립니다.

Hotplate Stirrer | AHS2001



제품 규격 및 사양

Model		AHS2001
Plate Dismension		180 x 180 mm
Stirrer Speed Range		80 to 1,500 rpm
Speed control resolution		5 rpm
Maximum Stirring Volume (H ₂ O)		20 Lit
Heating temperature		Max. 400 °C
Temp Accuracy		±0.3 °C
Temp Resolution	Display	0.1 °C
	Control	0.5 °C
Controller		Digital Fuzzy Controller with Touch Button & Rotary Switch
Timer		99hr 59min
Display		Full Digital LCD Display with Back-Light Function
Heating Output		600 W
Motor Type		Brushless Shade Motor
Permissible Temp. & Humidity		5°C ~ 40°C, <80%
Maximum Magnetic Stirrer Bar		50mm
Safety Devices		Over Temp Protection
Material	Plate	Ceramic Coated Aluminum
	Body	Powder Coated Aluminum
Power Consumption		650 W
Power Supply		1 Phase, AC 120V, 60Hz or AC 230V, 50/60 Hz

* 제품의 규격 및 사양은 사전공지 없이 변경될 수 있습니다. 발주전 문의 부탁드립니다.